



# Ceiling mounted

LED SERIES

LED EXAMINATION LIGHT | Focusmatic technology™

The DIASYS is an innovative LED examination light manufactured by SURGIRIS, made of aluminium. The DIASYS combines all the advantages of innovative LED technology with the most efficient energy management.



■ ON/OFF control intergrated to the light head.




# Technical data

LED SERIES

Reference	Diasys
Illuminated field diameter min/max at 1m	20 cm
Max. light intensity at 1m	65.000 Lux
Constant focal depth (L1+L2)	80 cm
Color temperature in °Kelvin	4500°K
Color rendering index CRI (Ra)	95
LEDs life time	>40 000 h
Feeding unit	100 to 240 V
Frequency	50 / 60 Hz

EN 60601-2-41

## Drawing



**Surgiris**  
taking care of you

SURGIRIS 80, rue de la Gare 59170 Croix - France  
Tel: +33 3 20 16 07 93 - Fax: +33 3 20 95 40 19  
contact@surgiris.com - www.surgiris.com  
Certified ISO9001- ISO13485

MADE IN FRANCE

