

C80 ICU Monitor




- CardioTec™ 12-lead ECG, Masimo SpO₂, AcuTec™ NIBP, RESP, TEMP, 12.1" LED touch screen, handwriting pen, EtCO₂ socket, Li-ion battery.
- Option: Dual IBP, EtCO₂, AG, Qcon, ICG, C.O., thermal printer, pediatric/Neonatal accessories.
- IPX1 level waterproof design.
- I-KLOK® intelligent alarm management for all parameters, auto-identification of alarm level.

C60 NICU Monitor



- Considering the feature of neonatal HR, we adopt the latest ExNeo® technology to guarantee the accuracy of measurement.
- World leading SpO₂ technology, Masimo SpO₂.
- Considering the low blood and low perfusion feature of neonate, we adopt Adap-DSP® system to provide more accurate measurement.
- Unique neonate measurement mode, soft & hard over-pressure protection.
- Bio-compatibility approved accessories for neonate.

COMEN



C50

Multi-parameter Patient Monitor

Shenzhen Comen Medical Instruments Co., Ltd.

Add: Floor 7, Block 5, 4th Industrial Park of Nanyou, Nanshan District, Shenzhen 518052, China

Tel: +86-755-2640 8879 2641 9446 Fax: +86-755-2643 1232

Website: en.comen.com E-mail: info@szcomen.com

FDA CE 0434



C50 Multi-parameter Patient Monitor Features



- 10.4" LED Touch Screen.
- Operating methods: Trim knob/ button, touch screen.
- Advanced patient file management system.
- High capacity storage, power-off data protection.
- Inclined screen for easy observation.
- Standard Configuration:**
5-lead ECG, Comen SpO₂, PR, NIBP, RESP, TEMP
- Optional Configuration:**
Masimo/ Nellcor SpO₂, Dual IBP, EtCO₂, Thermal Printer, C.O., ICG, AG



- Mounting solutions: vertical hanging, wall mounting, rolling stand.



- Complied with FDA, CE standards, integrated modular design and aluminum alloy structure guarantee product's safety & stability.
- High standard design, splendid appearance, best crafts work.



- IPX1 waterproof design, good performance in various environments.



- Fixed handle, easy to carry



- 4000mAh Li-ion battery, above 4 hours working time



- Brand-new interface design



- Handwriting pen, USB, network port, VGA interface.

- I-KLOK® Intelligent alarm management, automatic identification of alarm level. Self-adjust proper alarm time to reduce false alarms.

- ARR analysis, RESP Apnea Alarm, OxyCRG.

- CMRR≥105 dB, outstanding ECG anti-interference capability.

