

COMEN

Share With The World.



STAR5000

Fetal Monitor

Share the care brought by scientific and technical medical treatment when life starts

CE 0434



STAR5000 Fetal Monitor



Simple
operation mode



Superior
products performance



Excellent
system function



Powerful storage
and replaying system



Alarm system



Networking system



External printer

Standard

FHR / TOCO / FM / Inside Battery / 112mm thermal printer

Optional

External Printer / Twins / Waterproof sensor

Features

- **Simple operation mode**
Simple mode design, 12.1" color TFT touch screen, easy for observation
- **Superior products performance**
Advanced ARM system, convenient for maintenance and upgrade
- **Excellent system function**
Adjustment of FHR, TOCO waveform, parameter color, FM measurement auto/manual, automatic fetal status scoring system, more clinical information for better decision
- **Powerful storage and replaying system**
120 hours all-parameter displaying, 1 hour voice storage
- **Alarm system**
Intelligent 3-level audio/visual alarm, more accurate analysis
- **Networking system**
USB socket for software upgrade, ethernet socket supporting wire/wireless central monitoring for maximum 128 units
- **External printer**
Internal & external printer are available, improving quality & reducing cost



● 112mm wide thermal printer, paper rolls, Z-fold paper are available

● Hidden handhold, convenient to carry

● 12.1 inch LCD touch screen, convenient and efficient

● External printer is available



- Portable wallmount solution
- 12- crystal pulse wave sensor, for capturing more accurate FHR signal
- Main menu interface
- Patient information input and scoring interface



STAR5000D Fetal Monitor

- 6.5' color TFT touch screen, convenient for observation & movement.
- 152mm wide-line printer, excellent printing view.
- Automatic detection transducers, easy and flexible connection.
- 12 crystal wide pulse wave waterproof transducer, suitable for water birth.
- Advanced patient files management system, support handwriting input.
- TempPom timing monitoring function, avoid over time monitoring.
- Intelligent alarm management, automatic-identification of alarm level, self-adjust proper alarm time to reduce false alarms.



STAR5000E Fetal Monitor

- 5.7" folding-up LCD screen, easy for observation.
- 112mm printer, easy for paper replacement and observation.
- Tempom timing monitoring function, avoid overtime monitoring.
- Manual/automatic FM counting, automatic fetal status scoring system.
- Advanced patient files management system, support handwriting input.

Shenzhen Comen Medical Instruments Co., Ltd.

Add: Floor 7, Block 5, 4th Industrial Park of Nanyou, Nanshan District, Shenzhen City 518052, China

Tel: +86-755-2640 8879 2641 9446 Fax: +86-755-2643 1232

Website: www.szcomen.com E-mail: info@szcomen.com