

# EyeMax 2

*“Optimum comfort & protection  
during phototherapy treatments”*




## ▶ EyeMax 2 Overview

Our **EyeMax2** was designed based on valuable feedback from our customers. We have improved our design to assure optimum comfort and protection for our patients during phototherapy treatment. The single-piece, wrap-around design attaches at two points on either side of the baby's head and can be independently adjusted to prevent unwanted movement and allow a perfect fit. Our special UV light eye protection pad has been ergonomically designed to prevent any "light leaks

### Part Numbers

Regular (Comes in pack of 20) . . . . . R300P01  
 Premie (Comes in pack of 20) . . . . . R300P02  
 Micro (Comes in pack of 20) . . . . . R300P03




**1. Velcro strapping**

*Adjustable band for optimum comfort a perfect fits.*

**2. Eye Pads**

*Bigger eye pads for less pressure and better UV protection.*

**3. Latex Free** 

*Free from irritant contact.*

---

**Specifications**

Regular:.....	12.6" - 14.9" (Occipital Frontal Circumference)
Preemie:.....	10.4" - 12.6" (Occipital Frontal Circumference)
Micro:.....	7.87" - 10.4" (Occipital Frontal Circumference)